



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## SHREWSBURY SM3370 Ultrasound therapy device (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in UK,
- Power supply: 230 V,
- Emergency power supply: through a new built-in 12 V 2 Ah battery,
- Frequency: 50/60 Hz,
- Power: 40 VA,
- Technical specifications:
  - Dual frequency: 1 (continuous operation) / 3 (pulse operation 1:1, 1:3, 1:9) MHz,
  - Continuous output power range: 0 - 1.5 W/cm<sup>2</sup>,
  - Dual frequency head: 5 cm<sup>2</sup>,
  - Therapy timer range: 0 - 20 min,
  - 3 pulse options,
  - Built-in head holder,
  - Automatic power reduction without load,
  - Can be connected to combination therapy,
- Included:
  - Power cord,
  - Carrying bag,
  - Head,
- Dimensions:
  - Devices: 27 x 28 x 8 cm,
- Weight: 2.4 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

