



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DRAGER BABYTHERM 8004 Neonatal incubator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- DRAGER Babytherm 8004 is a neonatal intensive care station, designed for heat therapy for both babies born at term and premature babies,
- The device can be used in operating rooms, neonatal wards, children's wards, maternity wards and pediatric intensive care units. The device can be operated by clinical staff or nursing staff under medical supervision,
- Babytherm's mission is to provide the neonate with an optimal thermo-neutral environment, that is, to prevent thermal stress resulting from overheating or hypothermia,
- Ceramic heating elements and special gold reflectors ensure even heating of the entire surface of the mattress, so the newborn will always be properly warmed, regardless of the position on the mattress,
- The ability to simultaneously observe central and peripheral temperatures provides transparent information about the stability of the patient's temperature, which facilitates rapid detection of thermal stress,
- Integration of information from all medical devices at the care station ensures the highest optimization of work by minimizing errors, and saving time during basic procedures,
- The SmartSWIVEL mechanism ensures that the patient is warmed evenly at all times, even when the radiant heater is moved out of the way to widen the field of access to the patient. As a result, the baby is not subjected to the stress of hypothermia during, for example, X-rays or other procedures requiring free space directly above the crib,
- Two sensors ensure precise measurement of skin temperature also during kangarooing,
- The working light is equipped with day and night modes. Two brightness levels ensure proper illumination of the crib at any time,
- Incubator equipped with:
  - Control unit:
    - Located directly under the radiant heater, it is accessible from any side of the incubator, making it extremely easy and convenient to change settings,
    - Information is displayed on a large, bright digital display, visible from all sides,
    - Understandable alarm messages displayed on the control panel allow the operator to quickly remove the cause of the alarm,
  - Radiator with lateral rotation along with automatic tilt,
  - Precise DRAGER SCALE located under the control unit with a permissible patient weight of up to 8 kg,
  - An indicator light that allows the operator to quickly locate the alarming incubator,
  - Electric height adjustment of the incubator with foot switches within 30 cm,
  - Two available rotating cabinets with dimensions: 32 x 32 x 12 cm along with a rotating shelf with a maximum load capacity. 3 kg,
  - Hanger with 3 fluid container holders including height adjustment of 48 cm,
  - 6x airlock allowing wiring of sensors and hoses to the inside of the incubator,
  - One-sided tilting top cover allowing direct contact with the newborn,
  - Lever-operated mechanism that allows for Trendelenburg positioning on both sides within 20°,
- Mobile on 4 casters with built-in brake,
- Body temperature sensors - not included (can be retrofitted at additional cost),
- Selectable menu in language: EN, DE,
- Included:
  - User manual in Polish (PDF),
  - 2 x power cord,
  - Reusable mattress with medical approval,
- Dimensions:
  - Overall: 130 x 70 x 190/220 cm,
  - Mattress: 72 x 48 cm,
- Weight: 135 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

