



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VERMEIREN EAGLE 625 Patient transport elevator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2018),
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V 2.9 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 230 V,
- Lifters for easy transport and lifting of patients. Model 625 designed for patients weighing up to 175 kg . The lifts operate automatically by remote control , equipped with high-quality LINAK actuators guarantee smooth and trouble-free operation,
- Features:
 - The design allows easy and quick folding for transport without using any tools,
 - Models equipped with locking casters,
 - The running frame is movable and allows infinitely adjustable spacing for maximum stability,
 - Lifters equipped with an emergency button,
 - Headrest suspension made of durable water permeable material - quick drying,
 - The design of the lifts allows lifting patients from the bed, from the wheelchair,
- Patient lift equipment allows it to be used in buildings for patients affected by, for example:
 - Paralysis,
 - Loss of limbs (leg amputation),
 - limb damage / deformity,
 - stiff or damaged joints,
 - Heart and circulatory failure,
 - elderly people,
- Technical specifications:
 - Lowest arm position: 76 cm,
 - Highest arm position: 184 cm,
 - Minimum distance between the legs of the base: 68 cm,
 - Maximum distance between base legs: 122.5 cm,
 - Number of operations: 40 ~ battery charging,
- Included:
 - User manual in Polish(PDF),
 - Power cord,
 - 3 COMFORT carriers (Up to 125 kg):
 - 2x polyester - size L
 - 1x polyester mesh - size L,
- Dimensions:
 - Base length: 129 cm,
 - Base height: 11 cm
- Maximum load capacity: 175 kg,
- Weight: 42 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

