



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## ASPEL ASCARD B56 EKG Unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Emergency power supply: through 1 built-in 12 V 0.8 Ah battery -possible to replace with a new one for an additional fee,
- Technical specifications:
  - Frequency: 50/60 Hz,
  - ECG signal: 10 standard leads,
  - Real-time 3, 6, 12 - channel printing (paper width: 11.2 cm),
  - Test sensitivity, 2.5/ 5/ 10/ 20 mm/mV,
  - Built-in 2 x 24 character LCD screen,
  - Foiled, touch and readable control panel,
  - Signaling of bad connection of individual electrodes,
- Selectable menus in the following languages: PL, ENG, GER , FR, IT, ESP, HR, TUR, CS, RUM, SK, LT, SCC,
- Included:
  - User manual in Polish/English (PDF),
  - Power cord,
  - ECG 10-lead KEKG cable,
  - EKG clip limb electrodes - 4 pcs,
  - EPP suction cup electrodes - 6 pcs,
- Overall dimensions:
  - 25 x 19 x 6 cm,
- Weight: 2 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

