



## **OLYMPUS LUCERA CLV-260SL Light source (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 500 VA,
- Bulb: Xenon MD-631 300 W,
- Technical specifications:
  - 17-stage light intensity adjustment,
  - Lamp life: 500 h. of continuous use,
  - Foiled control panel,
  - 4 levels of pressure switching (off, low, medium, high),
  - Frequency fluctuation: +1/- 1Hz,
  - Diaphragm pump,
- Included:
  - Power cord,
- Dimensions: 38 x 53 x 16 cm,
- Weight: 16 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
  - The warranty period does not include bulbs and other light sources installed in the device,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

