



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER OXYLOG 1000 Ambulance ventilator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- The Oxylog 1000 ventilator is designed to ventilate patients in emergencies. Designed for use in the field, the Oxylog 1000 is one of the most multifaceted devices of its kind, distinguished by its intuitive operation, robustness and mobility,
- The respirator has visual and audible alarms to control both airway pressure and supply pressure, contributing to patient safety,
- The power supply for all functions of the device is pneumatic, so the respirator is independent of the power source. Thanks to its design, the Oxylog 1000 features outstanding reliability, which is so important to face the harsh conditions that accompany rescue ventilation,
- Technical specifications:
 - Built-in pressure gauge: from -10 to 80 mbar,
 - PV indicator - supply pressure,
 - Ventilation frequency (Freq): 5 - 35 1/min,
 - Minute volume (MV): 3 - 20 l/min,
 - Maximum airway pressure (Pmax): 25 - 55 mbar,
 - Mixing switch: 100% O2 or 60% O2,
- Includes:
 - User manual in Polish (PDF),
 - Patient breathing hose with valve,
 - Oxygen supply hose,
- Dimensions: 22 x 21 x 12.5 cm,
- Weight: 3.3 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



PRAYSTON