



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



CHATTANOOGA INTELECT MOBILE 2776 Ultrasound therapy device (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Mexico,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 75 VA,
- Main product features:
 - User protocols allow the operator to set, save and change the parameters of each program to suit specific patient needs,
 - Memory for up to 10 user programs,
 - Lightweight and innovative design gives you the ability to place the unit on a cart, countertop, wall or use it as a mobile device,
 - Electronic Signature allows the system to be automatically calibrated to any head size,
 - Efficient approach to treatment setup using pre-set parameters,
 - The ultrasound therapy device features:
 - 5" backlit LCD display with adjustable brightness for better visibility in low-light conditions,
 - Ergonomic head with heating function and contoured grip at 20° angle,
- Specifications:
 - Operation mode:
 - Pulsed (10, 20, 50%),
 - Continuous (100%),
 - Frequency: 1 or 3 MHz,
 - Pulse frequency: 16, 48, 100 Hz,
 - Pulse duration: 1, 2, 5 ms ± 20%,
 - Amplitude:
 - 0 - 2.5 W/cm² in continuous modes,
 - 0 - 3 W/cm² in pulsed modes,
 - Output:
 - Pulsed: 1 MHz signal modulated 100% by 100 Hz rectangular wave with selected duty cycle,
 - Continuous: 1 MHz nominal signal active as long as timer runs,
 - Maximum treatment time: 30 min,
 - 5 cm² head model 27335:
 - Frequency: 1.0, 3.3 MHz,
 - Power: 0 - 4 W,
 - Effective radiation area (ERA): 4.0 cm²,
 - Beam type: collimated,
 - Maximum beam non-uniformity ratio (BNR): 5:1,
 - Head suitable for use in water environment (IPX7),
 - Coupling: 1 - 4 W,
- Selectable menus in EN, TR, NR, SE, NL, NE, PT, IT, DE, FR, ES,
- Included:
 - User manual in English (PDF),
 - Power cord,
 - Head holder,
 - Head 27335,
- Dimensions:
 - Camera: 26 x 16 x 28 cm,
 - Head:
 - Length: 23 cm,
 - Cable length: 161 cm,
- Weight: 2.2 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

