



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **HILL-ROM ISOLETTE C2000 Neonatal incubator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- HILL-ROM Isolette C2000 is a neonatal incubator that provides a high level of thermoregulation and an ideal micro-environment for the development of newborns,
- Equipped with a cream-colored medically approved mattress,
- Has built-in systems:
  - Maintaining constant temperature,
  - Regulation of air humidity,
  - Alarm,
- Equipped with:
  - 6 nursing openings for easy access to the patient,
  - Distilled water container,
  - Rails for attaching additional accessories,
  - Grommet seal for wires and drains - 10 pieces. Socket for connecting two skin temperature sensors (one temperature sensor included),
  - Built-in patient weight (accuracy in 10 gram increments) - measurement up to a maximum of 7 kg,
  - Built-in control panel screen,
- Technical specifications:
  - Infinite manual change of the recliner angle to Trendelenburg and anti-Trendelenburg positions within 15°,
  - Air temperature: 20 - 37° C,
  - Skin temperature: 34 - 37 ° C,
  - Pressure: 40 - 150 psi,
  - Heating time: < 35 min,
  - Fresh air flow rate: 30 L / min,
  - Oxygen: 21 - 65%,
  - Noise level inside the chamber: < 47 dB,
- Selectable menu in language: PT, SE, FN, GR, PL , NR, DA , NL, JP, IT, ES, DE, FR , EN,
- Mobile on 4 casters Ø 12 cm all with brakes,
- Includes:
  - User manual in English (PDF),
  - Power cord,
- Dimensions:
  - Mattress: 37 x 72 x 2 cm,
  - Chamber: 45.5 x 83.5 x 49 cm,
  - Overall: 70 x 105 x 140 cm,
- Weight: 79 kg,
- Cat 04,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

