



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EDAN MT-803 Medical cart (Ex-Display)

- New (ex-display),
- Technical condition: perfect,
- Visual condition: perfect,
- Real photos of the product,
- Made in the UK, (2021),
- Medical cart MT-803 EDAN on a mobile base with 4 braking wheels,
- Built-in handle for free movement of the cart,
- Pull-out top and bottom drawer for accessories,
- Dimensions:
 - Top shelf: 44.5 x 45 cm,
 - Lower shelf: 37 x 31 cm,
 - Pull-out top shelf: 22 x 25 cm,
 - Pull-out bottom shelf: 29 x 31 cm,
 - Total: 44.5 x 50 x 92.5 cm,
- Weight: 22.5 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

