



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## OPTA KATOWICE Trial lens set (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Set of lenses of different strengths, which, together with the test frame, are used for eye examination and selection of appropriate glasses,
- Technical specifications:
  - 26 pairs of spherical lenses with positive and negative power:
    - 0.25 D - 3.00 D, in 0.25 D increments,
    - 3.50 D - 5.00 D, in 0.50 D increments,
    - 6.00 D - 10.00 D, in increments of 1.00 D,
    - 10.00 D - 12.00 D, in increments of 2.00 D,
    - 13.00 D - 14.00 D, in 1.00 D increments,
    - 16.00 D - 20.00 D, in 4.00 D increments,
  - 12 pairs of positive and negative power spherical lenses:
    - 0.50 D - 2.00 D, with plug in 0.25 D increments,
    - 2.50 D - 5.00 D, in 0.50 D increments,
- Included:
  - Trial frame:
    - Has spaces for 3 pairs of lenses,
    - Diameter of matching lens: Ø 38 mm,
  - Spacing adjustment,
- Missing 4 filters,
- Overall dimensions: 50 x 29 x 7.5 cm,
- Weight: 2.5 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

