



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ENVITEC MYSING O Oximeter (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2020),
- Power supply: 230 V,
- Emergency power supply: via 1 built-in Li-ion 3.6 V 2900 mAh battery (can be exchanged for a new one at additional cost),
- Frequency: 50/60 Hz,
- Power: 5 V,
- Technical specifications:
 - TFT color display:
 - 240 x 320 px,
 - 2.8' inch,
 - Measurement range: 0 - 100% oxygen,
 - Display accuracy: 0.1% oxygen,
 - Accuracy: <1% of O2 volume / 8 hours,
 - Sensor type: OOM111,
 - Operating time on one charge: 24 h,
 - IP 54 immunity certification (dust and water resistant),
 - IK 05 class impact resistance,
 - Software version: 01.95.00-WE.03,
- Selectable menu in the following languages: DE, EN, FR, ES, IT, NL, PL,
- Included:
 - User manual in English (PDF),
 - Power cord,
 - New oxygen sensor,
 - Adapter,
- Dimensions: 1.6 x 7.2 x 3.9 cm,
- Weight: 330 g,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

