



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDIMA S1 Infusion (syringe) pump (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Emergency power supply: through 1 new built-in battery: 9.6 V / 2 Ah,
- Frequency: 50/60 Hz,
- Power: 20 VA,
- Medima S1 is a single-channel syringe infusion pump for precise fluid dosing,
- The pump features:
 - Automatic syringe recognition and grasping,
 - Intuitive operation,
 - Visual and audible alarms with adjustable volume and beep type,
 - Backlit, easy-to-read graphic display,
 - Menu in language: Polish,
 - Foiled control panel in language: English,
 - Tripod mounting bracket,
 - Wide selection of units,
 - Night mode; reduced display brightness; programmable sound volume, programmable start and end time,
 - Simple connection and disconnection from docking station (without interrupting infusion),
- Infusion types:
 - Continuous,
 - Intermittent,
 - Profile,
- Syringes supported:
 - 5 mL - B. BRAUN OMNIFIX, BD PLASTIPAK, CODAN, ERG Kłobuck, KD-JECT III, NIPRO EUROPE, ONCE, TERUMO,
 - 10 mL - B. BRAUN OMNIFIX, BD PLASTIPAK, CODAN, ERG Kłobuck, FRES INJECTOMAT, KD-JECT III, MONOJECT, NIPRO EUROPE, ONCE, TERUMO,
 - 20 mL - B. BRAUN OMNIFIX, B. BRAUN PERFUSOR, BD PLASTIPAK, CODAN, ERG KŁOBUCK, KD-JECT III, MONOJECT, NIPRO EUROPE, ONCE, POLFA LUBLIN, TERUMO,
 - 30 mL - B. BRAUN OMNIFIX, BD PLASTIPAK, CODAN, MONOJECT, NIPRO EUROPE, ONCE, TERUMO,
 - 50 mL - B. BRAUN OMNIFIX, B. BRAUN PERFUSOR, BD
 - 50 mL - B. BRAUN OMNIFIX, B. BRAUN PERFUSOR, BD PERFUSION, BD PLASTIPAK, CODAN, DIPRIVAN, ERG KŁOBUCK, FRES INJECTOMAT, JANPOL, JIANGSU, KD-JECT III, MARGOMED, MONOJECT, NIPRO EUROPE, ONCE, POLFA LUBLIN, SHANCHUAN, TERUMO,
- Specifications:
 - Infusion rate: 0.1 - 2000 mL / h,
 - Infusion volume: 0.1 - 10000 mL,
 - Infusion time: 1 s - 99 h,
 - Bolus dose: 0.1 - 50 mL / h,
 - Bolus infusion: 2 - 2000 mL / h,
 - KVO: 0 - 5 mL / h,
- Included:
 - Power cord,
 - User manual in language: Polish (PDF),
- Dimensions: 36.5 x 23 x 12 cm,
- Weight: 2 kg,
- Cat 04,
- The device has a display defect in the form of several longitudinal stripes on the matrix of the screen display - the shortcoming does not affect the comfort of use,
- It has an up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

