



## **HAAG-STREIT BERN Ophthalmic table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Switzerland,
- Ophthalmic table with manual height adjustment is an elegantly and solidly made table for eye examination equipment such as slit lamp, auto-refractometer, tonometer, perimeter, etc.
- Technical specifications:
  - Mobile table on 4 casters (2 braking),
  - Manual adjustment in the range: 76 - 90 cm,
  - Pull-out drawer,
- Dimensions:
  - Overall: 54 x 35 x 76 / 90 cm,
  - Pull-out drawer: 13 x 31 cm,
- Weight: 41 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

