



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Maico Quik Tymp 2 Tympanometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Maico Quik Tymp 2 is designed for both tympanometry and stirrup reflex threshold measurement in adults and children,
- Technical Specifics:
 - Real-time graphical presentation of measurements,
 - Quick printout of measurement results,
 - Easy-to-read LCD screen with dimensions: 14 x 5 cm,
 - Built-in thermal printer with paper width: 11 cm,
- Menu in English,
- Keypad on the control panel in: English,
- Includes:
 - User manual in English(PDF),
 - Power cord,
 - Power supply,
 - Tympanometer head,
 - Earphones,
 - Patient button,
 - Carrying bag,
- Dimensions:
- Instrumentation bag: 40 x 30 x 20 cm,
- Weight: 6.5 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

