



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LEISEGANG 1H3 Colposcope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Leisegang 1H3 is a mobile colposcope for gynecological examination,
- Light source: halogen bulb 12 V / 20 W,
- Green filter,
- Mounted on a 3-part pantographic arm with a total length of about 130 cm,
- 3-segment arm adjustable in a large range of side-to-side rotation,
- Possibility to mount the colposcope arm on the wall and on a base / operator's table adapted to the arm,
- Infinite manual adjustment:
 - Height adjustment of the colposcope within 8 cm,
 - Light source intensity,
 - Pupil distance and eyepiece correction,
- The colposcope has traces of use that do not affect the comfort of use,
- Included:
 - Power cord,
 - Built-in light source,
- Dimensions: 60 x 130 cm
- Weight: 11.5 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

