



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARL ZEISS HUMPHREY 720 (735) Polometer (perimeter) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 400W,
- 12V/30W halogen bulb,
- Technical specifications:
 - Test type: static, single-eye and dual-eye driver tests,
 - Projector movement range: 90o,
 - Maximum pulse intensity: 10,000 ASB,
 - Pulse duration: 200 ms,
 - Canopy illumination: 31.5 ASB,
 - Test distance: 30 cm,
 - Features a built-in 12" monitor,
 - Hard disk and floppy disk drive,
 - Electric chin adjustment: up and down and side to side,
 - Lens holder,
 - Mobile on original Zeiss table with electric height adjustment from 69-100 cm,
 - Installed software from model 735,
- Selectable menus in EN, DE, ES, FR, IT, JP, PT, SE,
- Included:
 - User manual in English(PDF),
 - Power cords,
 - Keyboard,
 - Mouse,
 - LEXMARK 4513 printer,
 - Patient button,
- Dimensions:
 - Polometer: 56 x 46.5 x 58 cm,
 - Overall: 90 x 60 x 184 cm,
- Weight: 65 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

