



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BT-200 Fetal doppler (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Korea,
- Power supply via 2 built-in AA 1.5 V batteries,
- Frequency: 50/60 Hz,
- BT-200 is a portable and handy fetal heart rate detector. The device has adjustable volume for fetal heart rate, and a 3.5mm jack headphone jack. The LCD screen displays fetal heart rate and battery status. The ultrasound head operating at 2 MHz improves the reading quality of the device, and increases the chance of detecting the baby's heart rate in the early weeks of pregnancy,
- Operating frequency (MHz): 2;
- Waterproof certification: IPX7,
- Included:
 - User manual in English(PDF),
 - Original carrying case,
- Dimensions: 7.5 x 2.6 x 12.8 cm,
- Weight: 190 g,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

