



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHIMMELBUSCH Sterilization container (small) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Specialized cylindrical shaped container is used to sterilize dressing materials (compresses, gauze pads) in the autoclave to store them later as "sterile" materials for medical activities.
- For this purpose, the instrument has a circular shape, and on one side there is a perforation that allows steam to reach the sterilized materials. Moreover, the sterilization tool has a locking lid. The Standard sterilization can is made of high-quality surgical steel with very good anti-corrosion properties. It can be sterilized repeatedly.
- Dimensions:
 - Inner diameter: 15 cm,
 - Height: 15 cm,
- Weight: 1.5 kg
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

