



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **BQE DYNAMIC Medical chair (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the Netherlands,
- Dynamic is an ergonomic medical chair that activates and corrects the sitting position. It has a unique patented scissor mechanism between the seat and backrest. The lumbar support automatically moves backward and adjusts to the hollow contours. The natural S-shaped anatomical shape of the spine remains even when sitting passively. The small lumbar support provides optimal freedom of movement and activates spinal stretching,
- New medical approved upholstery in black,
- Perfect for doctors' offices and laboratories or, for example, as a dental chair,
- The base of the cross is made of chrome-plated steel, equipped with five double casters,
- Dimensions:
  - Seats: 37.5 x 43.5 x 9 cm,
  - Overall: 60 x 60 x 76 cm,
- Maximum load capacity: 110 KG,
- Weight: 11.5 KG,
- Cat 02,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

