



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **G.HUG GMBH ARTROMOT K CPM rehabilitation splint (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Artromot K is a device for continuous passive mobilization of the knee and hip joint, providing the possibility of improving patients who remain in bed,
- Allows individual selection of all parameters and a wide range of their adjustment,
- Functions:
  - Continuous passive motion therapy (CPM) of the knee and hip joint,
  - Knee joint flexion range: 0° / 110°,
  - Adjustable pause range: 0 - 20 sec,
  - Can be used for the left and right leg and fully adapted to the patient's anatomy (adjusting the device to the length of the femur and tibia),
  - Control by potentiometers (flexion, extension, speed, pause),
  - Setting of any flexion and extension angles, exercise speed, pause time during flexion or extension,
  - Possibility of active sole or dorsal foot flexion exercises,
  - Control panel in German,
- Included:
  - Power cable,
- Dimensions:
  - Rails: 100 x 37 x 60 cm,
  - Overall: 100 x 40 x 130 cm,
- Weight of set: 28 kg,
- Comes with a current Technical Passport issued,
- Warranty:
  - 6 months for domestic market - Poland,
  - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

