



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER CATO EDITION Anesthesia machine (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Has one port for evaporator with common Dräger DW-2000 connector systems, no evaporator included - can be retrofitted at additional cost,
- Supplied with gases from a central system: O2, AIR, N2O,
- Possibility to supply gases from additional O2 or N2O cylinders,
- Respirator driven by electric motor (electronically controlled),
- Built-in DRAGER gas measurement monitor with water trap:
 - With inspiratory and expiratory measurement of oxygen, nitrous oxide, carbon dioxide and 5 anesthetic gases: DES, ENF, HAL, ISO, SEV,
 - Audible and visual alarms,
 - Foiled control panel,
- Ventilation modes:
 - Manual, spontaneous , IPPV, PCV, SIMV,
 - I:E ratio 1:3 - 2:1, smooth,
 - PEEP 0:2 - 15 mbr,
 - Inspiratory volume 10 - 1400 ml,
 - Inspiratory rate 6 - 60/min. (IPPV),
 - 3 - 60/min (SIMV),
- Infinitely adjustable gas flow rate: O2, AIR, N2O, from 0.0 to 12 L / min, (digital displays and pressure gauges),
- Compact heated patient system,
- Built-in LED monitor with gas measurement function: O2, AIR, N2O, and control panel and ventilator function display,
- Equipped with a new O2 sensor,
- Easy to use and read display,
- Space for two gas cylinders,
- DRAGER INJECTOR SUCTION JAR:
 - DRAGER suction jar with a capacity of 700 ml suitable for phlegm extraction,
- Work surface with dimensions: 42.5 x 37 cm,
- Drawer with dimensions: 35 x 27 x 10.5 cm,
- Side rails for attaching additional accessories,
- Mobile, on casters, including 2 with brake,
- Can be retrofitted with a suction unit of your choice at an additional cost,
- Apparatus equipped with:
 - DRAGER INFINITY VISTA XL patient monitor/cardiac monitor,
 - Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Manufacture: German (2008),
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - Foiled and easy to read control panel,
 - Acoustic and visual alarms,
 - Liquid crystal TFT-LCD screen,
 - Screen resolution: 640 x 480 px,
 - Inputs:
 - Memory card,
 - Ethernet RJ45,
 - Sync,
 - RS-232,
 - Analog,
 - MultiMed,
 - NBP,
 - Wired modules:
 - ECG - Factory new 5 - lead ECG cable,
 - SPO2 - Factory new sensor + extension cable,
 - NIBP - Factory new cuff + air lead,
 - ECG - Factory new MULTIMED collection cable,
 - TEMP - No sensor included - can be retrofitted at additional cost,
 - Selectable menu in language: EN, DE,
 - Anesthesia machine equipped with a new O2 sensor
 - Included:
 - DRAGER VISTA XL monitor user manual in Polish(PDF),
 - Power cord,
 - Gas supply hoses with tips of customer's choice (AGA or DIN),
 - Patient breathing system,
 - Sampling line for gas measurement,
 - ECG - Factory-new 5-lead ECG cable,
 - SPO2 - Factory new sensor + extension cable,
 - NIBP - Factory new cuff + air line,
 - ECG - Factory new MULTIMED collection cable,
 - Dimensions: 96
 - Dimensions: 96 x 73 x 176 cm,
 - Total weight: 155 kg,
- In the domestic market (Poland), in the case of installation of the device within 30 days from the date of delivery of the equipment, installation and training are free, after exceeding this deadline these services will be charged,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

