



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAS SPECIAL MAGNETOTHERAPY (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- MAS-Special is a digital magnetic field therapy device.
- The device was developed by therapists and doctors specifically for home use,
- Technical specifications:
 - Pure sinusoidal frequency controlled by quartz,
 - Microprocessor control,
 - Frequency from 1 Hz to 9999 Hz (adjustable from 1 to 99 Hz in increments of one, from 100 to 9999 Hz adjustable in increments of ten),
 - High output power (50 watts at 10 Hz),
 - Power adjustable in 1% increments,
 - 60 preset and freely selectable programs,
 - Clear display for easy use,
 - Easy operation with 16-piece keypad,
 - Special lockable socket for connecting mat and cushion,
 - 10 coils anatomically arranged in the mat,
- Included:
 - Power cord,
 - Round applicator,
 - Applicator mat,
- Dimensions:
 - Mat: 45 x 65 x 15 cm,
 - Pillow: ø 20 cm,
 - Device: 26 x 26 x 12 cm,
- Weight: 15 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

