



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED ŻYWIEC SG 44.0 Gynecological segment with attachment and bowl (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Upholstery with medical certificate in black,
- Gynecological segment with attachment and bowl FAMED ŻYWIEC SG 44.0 is used to extend the seating segment in the table,
- It consists of a polyurethane mattress,
- Structure made of stainless steel with matte finish,
- Double-sided rails for attaching additional accessories,
- Urology attachment with segment attachment with 270° adjustable horizontal rotation,
- Bowl can be rotated 360°,
- Bowl with a capacity of 4 l,
- Included:
 - Urology attachment with and new waste bowl,
- Dimensions:
 - Mattress: 50 x 33 x 5 cm,
 - Bowl: 32.5 x 26.5 x 6.5 cm,
 - Total segment: 54 x 47 x 12 cm,
- Weight: 6 kg,
- Possibility to attach the segment to tables: SU-04, HYPERION,
- It has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

