



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VITALOGRAPH ALPHA 6000 Spirometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Ireland, (2012),
- Power supply: 230 V,
- Emergency power supply: via built-in new 7.2 V battery,
- Frequency: 50/60 Hz,
- Power: 30 VA,
- The Vitalograph ALPHA is an easy-to-use, lightweight desktop spirometer that provides the perfect combination of portability, speed and accuracy,
- Control panel keyboard in language: English,
- LCD screen with dimensions: 7.5 x 6 cm,
- Built-in printer for 11 cm wide paper,
- Flow detection principle: Fleisch type pneumotachograph,
- Back pressure: less than 0.1 kPa / l / sec. at 14 l / sec (compliant with ATS/ERS 2005),
- Volume detection: flow integration sampling at 100Hz,
- Maximum test duration: 90 seconds,
- Maximum displayed volume: 10 l ,
- Volume accuracy:
 - 3% or 0.05 l,
 - Maximum volume: 9.99 l,
 - Minimum volume: 0.01 l,
- Linearity: better than 3%,
- Accuracy at operating temperature range:
 - 10% or 0.3 l/s flow rate,
 - Maximum flow rate: 16 l / s,
- Minimum flow rate: 0.02 l / s,
- Operating temperature range:
 - ATS/ERS limits: 17-37 °C,
 - Design limits: 10-40 °C,
- Measures: VC, FVC, FEV1, PEF, FEF25-75% and more,
- On-screen visualization of volume/time and flow/volume graphs,
- Saves multiple initial tests for subsequent tests,
- Presents test information in a clear, easy-to-read format,
- Compact Flash memory card slot,
- Language in the control menu in: English,
- Included:
 - User manual in English (PDF),
 - Power cord,
 - Turbine with cable,
 - Original carrying case with dimensions: 46 x 37 x 11 cm,
 - 12 V power adapter,
 - USB cable,
- Dimensions: 30 x 25 x 8 cm,
- Weight: 3.7 KG,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

