



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TOPCON IT-1 Ophthalmic table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- TOPCON IT-1 ophthalmic table is an ophthalmic table with manual height adjustment by foot switch in the range; 67 - 93 cm,
- Automatic pushing of the table top to the highest position when there is no load,
- The table has pre-drilled holes in the table top prepared for mounting a slit lamp,
- Dimensions:
 - Tabletop: 55 x 37 cm,
 - Overall: 55 x 37 x 67/93 cm,
- Weight: 23 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

