



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



MEDELA THOPAZ PLUS Digital chest drainer (Ex-Display)

- New (ex-display),
- Technical condition: perfect,
- Visual condition: perfect,
- Real photos of the product,
- Made in Switzerland, (2022),
- Power supply: 230 V,
- Emergency power supply: via built-in 12 V Li-Ion battery,
- Frequency: 50/60 Hz,
- Thopaz sets new standards in chest drainage therapy. This compact system provides vacuum regulation very close to the patient's chest, based on the amount of air leakage. The device monitors air leakage and delivers the suction necessary to maintain negative pressure as determined and set by the attending physician. A digital display shows objective, real-time data and a 24-hour air leak graph, allowing easy monitoring of treatment progress. Setting up and operating the Thopaz is a breeze, and perfectly designed accessories further enhance the convenience of this lightweight system,
- Safety in the use of the Thopaz is ensured by alarms that alert hospital staff and call for quick intervention if an abnormality is detected. Thopaz continuously monitors and regulates negative pressure near the patient's chest, solving common problems experienced with unregulated drainage systems such as siphon formation, clogged or kinked drains. It is also possible to place the Thopaz at the patient's chest level.
- Advantages:
 - Unlimited patient mobility can promote faster recovery,
 - Objective real-time measurements and data showing changes in the size of the leak over time provide reliable information, allowing timely decisions on when to remove the drain,
 - Continuous patient monitoring provides peace of mind,
 - Easy setup and operation help hospital staff in their daily duties,
- Technical specifications:
 - Capacity: 5 L / min,
 - Vacuum:
 - 10 kPa,
 - 75 mmHg,
 - 100 mbar/cmH₂O,
 - Tolerance: +/- 15%,
- Included:
 - Power cord,
 - Docking station,
 - Fluid collection canister: 800 ml,
 - 4 pcs. - Patient drain,
 - Physical instruction manual,
- Dimensions: 25.5 x 9.5 x 23.3 cm,
- Weight: 1 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

