



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARL ZEISS ENT chair (Reconditioned)

- Reconditioned (used),
- Technical condition: perfect,
- Visual condition: perfect,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- New upholstery with medical certificate in black (upholstery can be exchanged for another one from our color chart for an additional fee),
- CARL ZEISS is a universal chair for medical examination,
- Technical specifications:
 - Electric height adjustment with foot control in the range: 44 - 60 cm,
 - Manual headrest adjustment in the range: 30 cm,
 - Possibility to manually adjust the headrest to the head position,
 - Seat rotates 360° around its own axis,
 - Patient side handles on both sides,
- Included:
 - Power cord,
 - Wired foot control,
- Dimensions:
 - Headrest: 21 x 19 x 6 cm,
 - Back rest: 42 x 58 x 6 cm,
 - Seat: 47 x 42 x 3.5 cm,
 - Overall: 60 x 60 x 122 cm,
- Weight: 62 KG,
- Cat 02,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

