



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RODENSTOCK Ophthalmic chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- New upholstery with medical certificate in light brown color (upholstery color can be changed, from our color chart for an additional fee),
- RODENSTOCK ophthalmic armchair sat on a sturdy base, together with the possibility of electric adjustment of the chair with a foot control,
- Technical specifications:
 - Electric height adjustment: 48 - 66 cm,
 - Manual adjustment:
 - Tilt of the chair in the range: 43°,
 - Height of the headrest in the range: 13 cm,
 - Headrest arm in range: 90°,
 - Head restraint: 180°,
 - Possibility of rotating the seat in the range: 360°,
 - Manually extendable center seat extension with lever,
 - Retractable patient step,
 - 2 retractable patient handles on both sides,
- Included:
 - Power cord,
 - Foot control,
- Dimensions:
 - Headrest: 21 x 17 cm,
 - Back rest: 48 x 50 cm,
 - Seat: 49 x 45 cm,
 - Step: 31 x 19 cm,
 - Overall: 66 x 77 x 128/146 cm,
- Weight: 85 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

