



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AGA URO-LIFT UUS-3060/EE Urology chair (New)

- New,
- Made in Germany,
- The chair consists of two parts,
- It has electric height adjustment, pelvic tilt, backrest adjustment by gas spring,
- Ergonomic upholstery of the back and seat provides a very high level of comfort. Integrating a foam wedge into the backrest eliminates the need for an additional headrest. There is a cutout in the seat section,
- Standard equipment:
 - A pair of standard rails attached to the side of the seat.
 - Two holes (right/left side) for a colposcope holder,
 - Grounded socket on the back of the column,
 - Cutout in the seat cushion,
 - Synchronous centrifugal sink (according to DIN EN 1717),
 - Paper roll holder screwed to the back,
 - Foot switches,
- Technical specifications:
 - Frame: Lacquered,
 - Seat dimensions:
 - Length: 1150 mm,
 - Back/seat section: 800 / 350 mm,
 - Width: 580 mm (640 mm with standard rails),
 - Height adjustment: 650 - 850 mm (electric motor)
 - Upholstery segments: 2,
 - Thickness of upholstery in the seat section: 85 mm,
 - Backrest upholstery thickness: 85 - 135 mm,
 - Seat section adjustment: +30° (electric motor),
 - Backrest adjustment: 0°/+54° (gas spring),
 - Trendelenburg position: 30°
 - Maximum load capacity of 200 kg.
- Has a current Technical Passport issued,
- Warranty 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

