



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AGA URO-MED UUS-3058/E Urology chair (New)

- New,
- Made in Germany,
- The chair consists of two parts,
- It has electric adjustment of pelvic tilt, backrest adjustment by gas spring,
- Armchair with a fixed height of 865 mm,
- Ergonomic upholstery of the backrest and seat provides a very high level of comfort. Integrating a foam wedge into the backrest eliminates the need for an additional headrest. There is a cutout in the seating area,
- The chair has been retrofitted with a folding footrest for easy chair entry,
- Standard equipment:
 - A pair of standard rails attached to the side of the seat.
 - Two holes (right/left side) for the colposcope holder,
 - Grounded socket on the back of the column,
 - Cutout in the seat cushion,
 - Synchronous centrifugal sink (according to DIN EN 1717),
 - Paper roll holder screwed to the back,
 - Foot switches,
 - Folding footrest,
- Technical specifications:
 - Frame: Lacquered,
 - Seat dimensions:
 - Length: 1150 mm,
 - Back/seat section: 800 / 350 mm,
 - Width: 580 mm (640 mm with standard rails),
 - Height adjustment: fixed height 865 mm,
 - Upholstery segments: 2,
 - Upholstery thickness in the seat section: 85 mm,
 - Thickness of backrest upholstery: 85 - 135 mm,
 - Seat section adjustment: +30° (electric motor),
 - Backrest adjustment: 0°/+54° (gas spring),
 - Trendelenburg position: 30°
 - Maximum load: 200 kg.
 - Folding footrest: 490 x 320 x 250 mm (W x D x H),
- Has a current Technical Passport issued,
- Warranty 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

