



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## TECH MED Rehabilitation couch (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Upholstery will be replaced with a new upholstery, in a color chosen by the customer, from our color chart,
- TECH MED rehabilitation couch is designed for rehabilitation treatments,
- It has a simple and stable frame made of powder-coated steel profiles,
- Headrest adjustable by gas spring,
- Paper hanger,
- Height adjustable with lower crank in the range: 47 - 96 cm,
- Headrest adjustment in the range: -70 / 0 / 45°,
- Included:
  - Paper hanger,
  - Headrest stopper,
  - Back segment stopper,
- Dimensions:
  - Overall: 66 x 92 x 47/96 cm,
  - Headrest: 66 x 46 cm,
  - Reclining segment: 66 x 144 cm,
- Weight: 60 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

