



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **CIM CART Monitor stand (Reconditioned)**

- Reconditioned (used),
- Technical condition: perfect,
- Visual condition: perfect,
- Real photos of the product,
- Made in Germany,
- CIM CART mobile stand with monitor mount,
- Technical specifications:
  - Built-in attachment handle for free movement of the stand,
  - Adjustable tilt of mounting arm within 25°,
  - VESA 75/100 standard,
  - Built-in tunnel to hide cables,
  - Mount compatible with monitors:
    - Drager: Infinity C500/700 | Explorer,
    - GE: Carescape B850,
    - Nihon Kohden: Life Scope G7/G9,
    - Spacelabs: Xprezzon display,
  - Stand placed on a mobile 5 - arm base with 2 braking casters,
- Dimensions: 57 x 134 cm,
- Maximum load capacity of the monitor: 7.5 KG,
- Weight: 12 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



