



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NOUVAG VACUSON 60 Surgical suction unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2015),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- The Vacuson 60 suction pump motor drives two pistons to create a vacuum. This achieves a high suction capacity of 60 liters per minute,
- The two-piston vacuum generation motor predestines the Vacuson 60 vacuum pump for extended applications,
- Specifications:
 - Precise and infinitely adjustable suction capacity,
 - Designed for continuous operation,
 - Low noise and vibration,
 - Robust, reliable and proven technology,
 - Wide Application,
 - Easy to use and maintain,
 - Automatic overfill protection (diverter valve),
 - Oil-free dual piston pump,
 - Carrying handles,
 - Possibility to connect 2 five-liter jars,
 - Possibility to connect an on/off pedal,
- Includes:
 - User manual in English (PDF),
 - Power cord,
 - Serres 1 L jar + cartridge,
 - Foot controller,
- Dimensions: 36 x 28 x 30 cm,
- Weight: 12 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

