



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **EMS MEDI-LINK 70A Short-wave diathermy (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the UK,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 200 VA,
- Memory backup battery: 3.6 V / 140 mAH,
- Technical specifications:
  - Emitting wave frequency: F= 27.12 MHz,
  - Peak power: 0 to 50 W,
  - Pulse repetition rate: 5 to 400 Hz,
  - Pulse width: from 20 to 400  $\mu$ s,
  - Treatment duration: 0 to 30 min,
  - Built-in display: 240 x 128 pixels,
  - Foiled touchscreen control panel,
  - Applicator placed on an articulated arm adjustable in multiple planes,
- The whole placed on a mobile 4-wheeled dedicated cart (2 braking wheels),
- Selectable menu in the following languages: ENG, FR, GER,
- Included:
  - User manual in English (PDF),
  - Power cord,
- Dimensions:
  - Overall: 70 x 45 x 105 cm,
  - Pull-out shelf: 23 x 26 cm,
- Weight: 22 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

