



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED-1 L13/62 Diagnostic and treatment light (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Bulb: 12 V / 35 W,
- Famed L13/62 is a single-focus diagnostic and treatment light,
- The light's canopy set on a spiral arm allows adjustment in many planes,
- The possibility of rotation in the range of 360°,
- It has a stable base,
- Infinite manual adjustment of the degree of focus of the light,
- The lamp gives shadowless light,
- Technical specifications:
 - Knob for intensity adjustment,
 - Manual on/off button,
 - Built-in power cord,
 - Height adjustment within 46 cm,
- Includes:
 - Power cord,
- Dimensions:
 - Canopy: Ø 20 x 18 cm,
 - Overall: Ø 30 x 157/203 cm,
- Weight: 14 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
 - The warranty period does not include bulbs and other light sources installed in the device,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

