



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **HAUS PAUSCH Massage chair (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz
- HAUS PAUSCH massage chair with the possibility of heating the back,
- Free massage intensity adjustment system with a knob,
- Built-in warm air blast system towards the back,
- Massage chair time adjustment,
- On/off switch for massage and back heating,
- Back angle adjustment with a switch in the range: 0 / 90°,
- Technical specifications:
  - Adjustable back massage in the range of 80 cm (operator can adjust the length section on which the massage is to move),
  - Ability to adjust the working time of the chair in the range of 15 minutes,
  - Tilt of the chair with switches in the range: 0 / 90°,
  - Manual adjustment of the inclination of the patient's seat within the range with the crank: 0 / 15°,
  - Free adjustment of footrest along and at an angle,
- Included:
  - Power cord,
  - Patient button,
  - Patient attachment with adjustable straps,
- Dimensions:
  - Overall: 60 x 127 x 155 cm,
- Weight: 165 KG,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

