



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BTL 5110 Laser therapy unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the UK,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 80 VA,
- Software: 7.07.0000R100,
- BTL-5110 is an advanced single-channel laser therapy device. The device is equipped with a computer with pre-installed programs, therapeutic encyclopedia and patient database, in addition, it is possible to self-program exactly the programs we need,
- Technical specifications:
 - Screen: 5.7" color touch screen,
 - Built-in patient database,
 - Number of independent channels: 1,
 - Number of simultaneously connected probes: 2,
 - Operation mode: Continuous and pulsed,
 - Dose: Infinitely adjustable 0.1 - 100 J/cm²,
 - Frequency: 0 - 10,000 Hz,
 - Therapy area: 0.1 - 100 cm²,
 - Fill factor: 35 - 90%, 100%,
- Application:
 - Increasing ATP levels,
 - Stimulation of protein synthesis including collagen,
 - Change in the structure and function of biological membranes,
 - Increased activity of certain enzymes and prostaglandins,
 - Modulation of immune processes,
 - Phagocytosis effect,
 - Increased cell motility and reproductivity,
- Indications:
 - Wounds and surgical scars,
 - Ulcers and transplant scars,
 - Diseases of the musculoskeletal system,
 - Various skin diseases,
 - Rheumatism and inflammatory conditions,
- Device placed on a dedicated mobile cart with a fixed shelf for accessories:
 - Cart on 4 casters (2 braking),
 - Selectable menu in ENG, FR, TUR, CHIN, CS, IT, KOR, GER, PT, NOR, PL, RUS, SWE, ESP, RUM, KZ,
- Included:
 - User manual in English (PDF),
 - Power cord,
 - 2x laser probe:
 - Red spot probe 50 mW / 685 nm,
 - Infrared spot probe 400 mW 830 nm,
 - 2x screen stylus,
- Dimensions:
 - Overall: 54 x 45 x 131 cm,
 - Unit: 23 x 21 x 40 cm,
- Total weight: 30 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility of extending the warranty up to 12-months for domestic and international markets for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

