



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ MEDI-MATIC 115.515 Gynecological chair with ANDROMEDA ELLIPSE Urodynamics (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Schmitz Medi-Matic is a chair used for gynecological and urological examinations and procedures,
- The chair is retrofitted with a urodynamics kit for examining the patient's urine,
- Medically approved upholstery in navy blue (upholstery can be replaced with new upholstery from our color chart at an additional cost),
- Smooth surface for optimal cleaning and hygiene,
- Infinite electric adjustment of the chair by foot switch:
 - Seat height in range: 74 - 95 cm,
 - Back segment angle in the range: - 16° to 43°,
- The chair has:
 - Paper space with width: 40 cm,
 - Rails for attaching additional accessories,
 - Manually adjustable foot rests,
 - Pull-out stainless steel bowl with dimensions: 29.5 x 32.5 x 6 cm,
 - Mountable leg rest with dimensions: 58 x 71 x 6 cm,
 - Mountable thigh rest on both sides,
- Equipped with:
 - ANDROMEDA ELLIPSE Urodynamics:
 - Membrane keyboard control (computer keyboard can be connected via PS jack),
 - Monochrome LCD screen with 320 x 200 px resolution,
 - Integrated VGA card for direct connection of a standard VGA computer monitor,
 - 200 dpi thermal printer - 110 mm print width,
 - The apparatus has the following measurement channels:
 - Pressure: 8 independent channels (Pves, Pura, Pabd, 5 defined pressures),
 - Flow: 2 independent channels (Qura1, Qura2),
 - EMG: 2 independent channels (EMG1, EMG2),
 - Infusion (filling) independent weight transducer,
 - Puller (puller),
 - Roller pump (for fluid administration only),
 - Ability to perform tests using urodynamics:
 - Urethral flow (automatic or forced start of the test),
 - Cystometry (filling phase) and Pressure/Flow (micturition phase),
- Includes:
 - SCHMITZ MEDI-MATIC chair:
 - Foot switch,
 - Power cord,
 - Waste bowl with dimensions: 32 x 30 x 6 cm,
 - Removable knee bolsters on both sides,
 - Side rails for attaching additional accessories,
 - Patient headrest,
 - Patient step,
 - ANDROMEDA ELLIPSE Urodynamics:
 - Power supply,
 - Micturition transducer,
 - Micturition funnel,
 - Swivel mount for micturition transducer,
 - EMG preamplifier,
 - IR remote control,
 - Cable for 1st pressure channel,
 - Cable for 2nd pressure channel,
 - Infusion transducer,
 - Camera mast for urogynecology chair,
 - Holder for three pressure transducers,
 - EURO clamp for attaching the handle to the micturition transducer,
- Dimensions:
 - Backrest: 58 x 87 x 5 cm,
 - Seat: 56 x 32 x 5 cm,
 - Base: 54 x 94 cm,
 - Overall: 122 x 90 x 133/154 cm,
- Maximum load capacity: 135 kg,
- Weight: 200 KG,
- Has a current Technical Passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 1 month for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

