



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARL ZEISS LR65844 Slit lamp (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Bulb: 12 V / 30 W,
- Frequency: 50/60 Hz,
- Technical specifications:
 - Eyepiece correction: (-5) / 0 / 5,
 - Magnifications: 5x / 8x / 12x / 20x / 32x,
 - Infinitely adjustable slit spacing, along with the ability to adjust the shape of the light beam with a knob:
 - Built-in filters:
 - Blue,
 - Green,
 - Infinitely adjustable:
 - Eyepiece spacing,
 - Patient's chin in range: 25.5 - 31.5 cm,
 - Side-to-side rotation,
 - Light intensity,
- Included:
 - Power supply with power cord,
 - Built-in fixation point,
- Dimensions: 32 x 40 x 63 cm,
- Weight: 12 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

