



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



REHA FUND Rehabilitation rotor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- REHA FUND rehabilitation rotor,
- The rotor is designed for rehabilitation exercises of the lower and upper extremities in traumatic-orthopedic and rheumatic diseases,
- Rotor made of high-grade steel and nickel alloys. Made of double reinforced frame,
- Has a wide range of resistance strength adjustment. Equipped with rubber straps to secure the foot,
- Dimensions: 42 x 50 x 30 cm,
- Weight: 2 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

