



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Forearm angle support (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Adjustable forearm support with the possibility of mounting on the mounting rail with a bone,
- Freely adjustable vertical tilt and side-to-side roll of the support.
- Infinite 360° horizontal adjustment,
- Infinite vertical adjustment with bone within 360°,
- Height adjustment within: 42 cm,
- Includes:
 - O.R. Direct mounting bone,
- Dimensions:
 - Backrests: 19 x 21 x 5.5 cm,
 - Overall: 19 x 21 x 62 cm,
- Weight: 2.5 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

