



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



RANGER 245 Blood and fluid warming device (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Capable of heating liquids at KVO flows up to 30,000 ml/h with special single-use cartridges,
- Heating power adjusted according to flow rate and inlet temperature,
- High heat conduction coefficient of the heating plates, between which there is a cassette with the flowing, heated fluid,
- Can be mounted on a stand with a base or on a regular drip stand. It is also possible to screw it to the operating table,
- Heating by means of dry heat,
- Specifications:
 - Fluid heating range: 33 - 41°,
 - Over-temperature alarm,
 - Alarm regarding insufficient temperature,
- Included:
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 19 x 25 x 11 cm,
- Weight: 3.4 KG,
- Has an up-to-date inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

