



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



IONOSON PHYSIOVAC 62 Electrotherapy unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- PHYSIOVAC 62 is a vacuum unit, controlled by a microprocessor. The device provides easy and convenient application of vacuum electrodes, especially in hard-to-reach areas,
- The orderly operating elements ensure extremely easy and comfortable work with the apparatus. The self-adjusting pump ensures quiet operation of the device and guarantees a constant value of the vacuum set by the therapist,
- Specifications:
 - Effective area: 2.5 cm² / 5 cm²,
 - Power output: 2.5 W / 5 W,
 - Average power density: 1 W/cm²,
 - Frequency: 1.0 MHz,
- Included:
 - Power cords,
 - Connection cable,
 - Vacuum electrodes:
 - 2 pcs. - 6 cm,
 - 2 pcs. - 8.5 cm,
 - Vacuum tubes:
 - 2 pcs. - Vacuum tube red,
 - 2 pcs. - Vacuum wire blue,
 - Vacuum connection blanking wire,
 - Ultrasonic head 2.5 cm²,
 - Operator button,
 - 4-outlet patient cable,
 - 6 pcs. - Plate electrode: 5.5 x 8.5 cm,
 - 1 pcs. - elastic bandage: 10 x 125 cm,
 - 2 pcs. - elastic bandage: 6 x 80 cm,
 - 7 pcs. - viscose cover for electrodes: 9 x 8 cm,
 - 2 pcs. - viscose pad for electrodes: ø 7 cm,
 - 2 pcs. - viscose pad for electrodes: ø 4 cm,
 - 2 pcs. - viscose covering: ø 11 cm,
 - 2 pcs. - viscose covering: ø 7 cm,
 - 1 pc. - point electrode tip: 3.5 cm,
 - 1 pcs. - point tip
 - 1 pc. - Point electrode tip: 2.5 cm,
 - Electrode holder (without buttons),
 - 2 pcs. - PHYSIOMED cover,
 - Electrotherapy film set,
- Device placed on a mobile dedicated cart (2 braking wheels),
- Dimensions:
 - Total: 55 x 45 x 80 cm,
- Total weight: 12 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

