



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HERATHERM THERMO SCIENTIFIC Incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2021),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1390 W,
- HERATHERM incubator with advanced protection has additional features that create the highest level of safety for your samples,
- Technical specifications:
 - Temperature range: 5°C TO 105°C,
 - Hysteresis: 0.5°C,
 - Capacity: 104 l,
 - Air circulation: double,
 - Alarms:
 - Audible door opening,
 - Failure to reach and over temperature,
 - Built-in RS-232 interface,
 - English language menu,
- Included:
 - English user manual (PDF),
 - 3 pcs tray with dimensions: 43 x 33 cm,
 - Power cord,
- Overall dimensions: 64 x 36 x 60 cm,
- Weight: 56 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

