



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS WM-NP1 Endoscope cart (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the UK,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Olympus WM-NP1 is a stable and sturdy cart suitable for carrying endoscopes and delicate medical equipment,
- It has:
 - Upper top with dimensions: 47 x 38 cm,
 - Lower top with dimensions: 47.5 x 56 cm,
 - 2 shelves with dimensions: 40 x 54 cm with manual height adjustment (every 1.5 - 2 cm) with Allen screws,
 - Handles for easy movement of the cart,
 - Endoscope holder with manual height adjustment,
 - Monitor holder with height adjustment within 20 cm (6 points every 3.5 cm),
 - Electrical installation, 12 power sockets for the medical equipment carried,
- Product features:
 - Power switch on the front of the cart,
 - Detachable rear panel of the cart for easy cable management,
 - Power cords aesthetically concealed in panels on the sides of the cart,
 - Corners of shelves protected with protective caps and edges over wheels with rubber,
- Mobile on four wheels with Ø 12 cm, two with brakes,
- Included:
 - 5 pcs. - power cord,
- Dimensions: 67 x 66 x 166/200 cm,
- Weight: 86.5 KG,
- Cat 03,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

