



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PRAISTON Hospital mattress (New)**

- New,
- Real photos of the product,
- Made in Poland,
- A medically approved mattress in gray color,
- It has zones of foam cut into the shape of waffles, allowing individual adaptation to the shape of the patient's body and reducing pressure and rubbing in the area of contact with the epidermis,
- Cut into the shape of waffles, the top foam facilitates air flow under the surface of the mattress, thereby counteracting the accumulation of moisture, which in turn helps reduce maceration and abrasion of the epidermis and dermis of the patient,
- It has a heel zone made of low-density foam, (soft) which helps to provide gentle support for the heels and ankles, making it easier for them to get off the mattress,
- It is transparent to x-rays,
- Does not contain latex,
- It is easy to clean and extremely hygienic - the cover can be machine washed or disinfected with a damp cloth,
- Fits most bed frames and adjusts to all patient positions in bed,
- Dimensions: 80 x 195 x 14 cm,
- Weight: 8 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

