



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TS-6830 Dental chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the PRC,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 800 VA,
- Medical approved upholstery in green color has traces of use (upholstery can be replaced with a new one from our color chart for an additional fee),
- Technical specifications:
 - Washable upholstery in green color,
 - Electric adjustment by wired remote control:
 - Height adjustment: 54 - 80 cm,
 - Back rest in the range of 55°,
 - Manual adjustment:
 - Headrest height: 16 cm,
 - Headrest tilt: 140°,
 - Armrest spacing: 180°,
- Included:
 - Power cord,
 - Remote control,
- Dimensions:
 - Backrests: 63 x 46 cm,
 - Backrest: 45 x 108 cm,
 - Headrest: 21 x 16 cm,
 - Overall: 62 x 170 x 68 cm,
- Maximum load capacity: 180 kg,
- Weight: 100 kg,
- It has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

