



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ANETIC AID QA4 Transport stretcher (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the UK,
- Power supply: 230 V,
- Emergency power supply: through 2 new built-in 12 V 4.5 Ah batteries,
- Frequency: 50/60 Hz,
- Anetic Aid QA4 is a functional transport stretcher that can serve as an operating table thanks to its numerous features,
- Upholstery with medical approval in black (upholstery can be replaced with a new one at an additional cost),
- Sleep surface made of double laminated reinforced HPL board,
- Mattress covered with waterproof antibacterial coating,
- Manual adjustment:
 - Headrest within $(-30)^{\circ} / 0^{\circ} / 25^{\circ}$,
 - Side tilt within 12° by crank,
- Electric control by remote control:
 - Height within 68 - 108 cm,
 - Backrest tilt angle in the range of 65° ,
 - Longitudinal movement in the range of 45 cm,
 - Changes in the patient's lying direction relative to the table column: Trendelenburg position within 20° and anti-Trendelenburg position within 12° ,
- Equipped with:
 - Handles on both sides of the headrest for convenient bed guidance,
 - Rails for attaching additional accessories,
 - A groove in the base for an oxygen cylinder,
 - 2 convenient brakes on both sides of the bed,
 - Directional wheel with lever,
 - 2 bumper disks,
 - Emergency handle that allows the user to immediately lower the backrest in case of an emergency,
 - Retractable railings on both sides,
- Table adapted for C-arm,
- Possibility to remove the leg segment and headrest,
- X - RAY translucent,
- Mobile, on 4 double wheels with $\varnothing 15$ cm,
- Includes:
 - User manual in English (PDF),
 - Power cord,
 - Controlling wired remote control,
 - Removable drip tray of 66 cm,
- Dimensions:
 - General: 211 x 80 x 68/108 cm,
 - Bed: 200 x 59 x 8.5 cm,
- Maximum load capacity:
 - Bed: 190 kg,
 - Head segment: 50 kg,
 - Leg segment: 25 kg,
- Weight: 153 KG,
- Cat 02,
- Has an up-to-date inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

