



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



IONTO-STERIL UL Hot air sterilizer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 690 W,
- IONTO-STERIL UL is a hot air sterilizer, is a fully safe and durable device. It meets all the requirements for compact devices,
- Safety is guaranteed by the function of automatic shutdown of the sterilization cycle including power off when the cycle is completed,
- The device can be used in veterinary and medical fields,
- Forced air circulation,
- Technical specifications:
 - Capacity: 4.3 L,
 - Infinitely adjustable temperature up to max. 180° C,
 - Has a timer switch adjustable from: 0 - 120 min,
 - Stainless steel housing, strong construction and long life of the device,
 - Dry sterilization,
- Includes:
 - Power cord,
 - Tray holder,
 - Tray - 1 piece,
- Dimensions:
 - Overall: 31 x 41.5 x 26 cm,
 - Tray: 18 x 28 cm,
- Total weight: 14 kg,
- Cat 02,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



