



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG VACUKLAV 24V Autoclave (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2500 W,
- Melag Vacuklav 24V is an automatic Class B autoclave equipped with MELAdem 40 water treatment station,
- Pre-vacuum and fractional vacuum produced by vacuum pump,
- Electronic LCD display,
- Device menu in German,
- Acoustic alarms,
- RS 232 printer socket,
- Has 6 operating programs:
 - Universal program: 134° C - metal tools, packaged,
 - Quick program: 134° C - metal tools, unwrapped,
 - Schon program: 121° C - packaged instruments, rubber, textiles, glass,
 - Disinfection: 105° C - unwrapped instruments,
 - Bowie & Dick test: 134° C, 2.2 bar 3.5 min - vapor permeation test to check the efficiency of the autoclave by the user,
 - Leakage test: vacuum test to check the leakage of the autoclave,
- Keypad on the control panel in the language: German,
- Included:
 - Factory new MELAdem 40 demineralizer, (water softening station):
 - Allows operator to quickly and cheaply obtain demineralized water. The device works on the principle of deionization and can be connected directly to the autoclave. The two-chamber system guarantees optimal consumption of the resin cartridge,
 - German manufacture,
 - Bed capacity: 210 l,
 - Capacity: 115 l/h at 700 µS/cm tap water,
 - Water purity: < 1-5 µS/cm²,
 - Dimensions: 32 x 35 x 15 cm,
 - 2 mixed resin cartridges for MELAdem 40,
 - Connection hoses - 2 pcs,
 - Perforated trays with dimensions: 18 x 37 x 1.5 cm - 2 pcs,
 - Power cord,
- Dimensions:
 - Chamber capacity: 19.5 L,
 - Chamber: Ø 22.5 cm, depth 45 cm,
 - Demineralizer: 32 x 35 x 15 cm,
 - Overall: 53 x 70 x 49 cm,
- Weight: 48 KG,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

